

● ABOUT US ●

Driven by the passion for innovation, we at Dr Trust endeavour to provide our customers with the latest medical inventions with an objective to promote good health and wellness all around the world. All the medical devices and health monitors provided by Dr Trust are supported by accurate, latest and ground breaking technologies, innovated at our headquarter in NY, USA. All our products adhere to the most stringent CE and FDA guidelines and are strongly recommended by doctors and health practitioners. Our products are designed in the utmost exemplary ways to ensure that their accuracy and convenience are unrivalled. The ease of their use and operation makes them even more suitable for users of all age groups.

Dr Trust strives to enhance the quality of lifestyle by providing with the most trusted and innovative health care and wellness products. Being a renowned global leader in health care products, Dr Trust ensures that our technically efficient team works dynamically and tirelessly to provide the best of the medical devices to our clients. The products that we have to offer are suitably designed for use at homes, laboratories and hospitals.

Our ground breaking solutions allow you to monitor your health in the easiest ways possible. In today's era when all of our lives are too hassled to handle, it becomes a bit difficult to pay attention to our health. But it has now become easier with the coming of the monitoring devices which can be conveniently used at homes and even on the go.

We bring to you a variety of best self medical devices, trusted and used by Doctors, medical professionals and home users all over the world.

● QUICK START GUIDE ●

Step 1

Slide the blood pressure cuff onto your upper arm and secure it so that it sits snugly about one inch above of your elbow.

Step 2

Simply push the START/STOP button and cuff begins to inflate with controlling deflation speed automatically.

Step 3

With the cuff inflation, look at the screen to get your blood pressure reading.

Step 4

Intended for giving very fast and reliable measurements, it gives result with Cuff Inflation.

Step 5

Hypertension Classification Indicator displays coloured results according to the range between which blood pressure values lie.

Step 6

Users can make settings for talking (ENGLISH/HINDI), memory and other features to make blood pressure monitoring convenient.



SMART SENSEI™ TECHNOLOGY

Dr Trust 360 App comes with our Smart Sensei™ technology that helps you take readings faster than ever. Take an image of your reading from your Dr Trust device and our app will automatically save it to your logs.

How does our Dr Trust 360 App work?



Download

Scan the QR Code to Download the Dr Trust 360 app and Register as a New User.



Click

Capture the results on your Dr Trust device screen through the camera built in the Dr Trust 360 app.



AI Driven Detection

The image is processed through our proprietary machine learning software, and converted into digits.



Monitor

Take note of your readings presented on meaningful graphs and monitor your health trends.




Share

Share your results seamlessly with family members, doctors and other health guardians.



drtrust360.com
The Science of Living



Note:

If you need to repeat the measurement, wait a few minutes before attempting a new reading.

Table of Contents

1. INTRODUCTION

- 1.1. Features of Dr Trust BP Smart
- 1.2. Important Information About Self-Measurement

2. IMPORTANT INFORMATION ON THE SUBJECT OF BLOOD-PRESSURE AND ITS MEASUREMENT

- 2.1. What causes high blood pressure?
- 2.2. Which values are normal?

3. THE VARIOUS COMPONENTS OF THE BLOOD-PRESSURE MONITOR

4. PUTTING THE BLOOD-PRESSURE MONITOR INTO OPERATION

- 4.1. Inserting the batteries
- 4.2. User Selection
- 4.3. Language selection, Time/Date settings



5. CARRYING OUT A MEASUREMENT

- 5.1. Before the measurement
- 5.2. Fitting the cuff
- 5.3. Start Taking BP Measurement
- 5.4. Discontinuing a measurement
- 5.5. Memory – storage and recall of the measurements
- 5.6. Memory full
- 5.7. Memory – cancellation of all measurements

6. APPEARANCE OF THE HEART ARRHYTHMIA INDICATOR FOR EARLY DETECTION

7. ERROR MESSAGES / MALFUNCTIONS

8. CARE AND MAINTENANCE, RECALIBRATION

9. BATTERY LIFE

10. SAFETY, CARE AND DISPOSAL

11. REFERENCE TO STANDARDS

12. TECHNICAL SPECIFICATIONS

13. REMARK

14. MANUFACTURER'S DECLARATION

15. CUSTOMER SUPPORT

16. ABOUT US

1. INTRODUCTION

1.1. Features of the Dr Trust BP Smart

The blood-pressure monitor Dr Trust BP Smart (with integrated time/date display) is a fully automatic, digital blood-pressure measuring device for use on the arm, which enables very fast and reliable measurement of the systolic and diastolic blood-pressure as well as the pulse frequency by way of the oscillometric method of measurement.

The device offers very high and clinical tested measurement accuracy and has been designed to provide maximum of user-friendliness.

Before using, please read this instruction manual carefully and keep it in a safe place. For further questions on the subject of blood-pressure and its measurement, please contact your doctor.

1.2. Important Information About Self-Measurement

- Do not use with neonatal patients.
- Do not intend to use with pregnant or pre-eclamptic patients.
- Keep product in proper packing. Any inadequate use will cause harmful injury to the patient or would affect the blood pressure due to connection tubing kinking.
- Too frequent measurements can cause injury to the patient due to blood flow interference.
- The application of the cuff over a wound can cause further injury.
- Avoid pressurization of cuff on any limb where intravascular access or therapy or an arteriovenous (A-V) shunt is present. It could result in injury to the patients because of temporary interference to blood flow.
- The patient is an intended user.



- Do not let the cuff put its pressure on the arm on the side of a mastectomy.
- Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring Medical Electrical (ME) equipment on the same limb.
- There is a need to check that operation of the automated sphygmomanometer does not result in prolonged impairment of patient blood circulation.
- Not intended to be used together with HF surgical equipment.

Electromagnetic Interference

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. IMPORTANT INFORMATION ON THE SUBJECT OF BLOOD -PRESSURE AND ITS MEASUREMENT

2.1. What Causes High Blood Pressure?

While the cause of high blood pressure in most people remains unclear, inactivity, poor diet, obesity, older age, and genetics-- can all contribute to the development of hypertension.

2.2. Which values are normal?

Blood pressure is too high, if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood-pressure is over 160 mmHg. In this case, please consult your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor.


Table for classifying blood-pressure values (unit: mmHg) according to World Health Organization

Range	Systolic Blood-pressure	Diastolic Blood-pressure	Measures
Blood pressure optimum	between 100 and 120	between 60 and 80	Self Check
Blood pressure normal	between 120 and 129	between 80 and 84	Self Check
Blood pressure slightly high	between 130 and 139	between 85 and 89	Consult your doctor
Blood pressure too high	between 140 and 159	between 90 and 99	Seek medical advice



Blood pressure far too high	between 160 and 179	between 100 and 109	Seek medical advice
Blood pressure dangerously high	Higher than 180	Higher than 110	Urgently seek medical advice!

3. THE VARIOUS COMPONENTS OF THE BLOOD-PRESSURE MONITOR



4. PUTTING THE BLOOD PRESSURE MONITOR INTO OPERATION

4.1. Inserting the batteries

- Insert the batteries (4 x size AA 1.5V), thereby observing the indicated polarity.
- If the battery power reduces less than 20%, battery warning icon appears.
- If the batteries are empty, battery warning icon appears indicating immediate replacement by new ones.




To Be Noted

- After the battery warning icon appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5V Batteries.
- Please remove the batteries from the device if the blood-pressure monitor is left unused for a long period.

4.2 User Selection

This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

- 1) Before measurement, make sure you set the user no. for the intended user. The unit can track results for 2 individuals (User 1, User 2) with their identity
- 2) Press the TIME button for at least 3 seconds. The display now indicates the set user. To confirm, press ON/OFF button.
- 3) Click the MEMORY button to select User 2.
- 4) We suggest the first person to take their pressure to be User 1.




4.3. Language Selection, Time / Date Settings

Language Selection

User can select either of the two languages i.e. English & Hindi, if talking required. Else user can also mute the unit.

- 1) Press & hold on/off button for 3 seconds. Default Language 1 (Language English) will appear.
- 2) To change language to Language 2 (Hindi), press memory button one time. Press this button twice to mute. To go back to English press setting button.




Setting the Time & Date

This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement can also be stored. After new batteries have been inserted, the clock originally displays DATE as 1-01 and TIME as 12:00. You must then re-enter the date and current time. For this, please proceed as follows.

1. Press the TIME button for at least 3 seconds firstly, user icon will blink. Then press TIME button again, the display now indicates the set year, during which the four characters blink.
2. The correct year can be entered by pressing the MEMORY button.
3. Press the TIME button again. The display now switches to the current date, during which the first character (month) blinks.
4. The corresponding month can now be entered by pressing the MEMORY button.



5. Press the TIME button again. The last two characters (day) are now blinking.
6. The corresponding day can now be entered by pressing the MEMORY button.
7. Press the TIME button again. The display now switches to the current time, during which the first character (Hour) blinks.
8. The corresponding hour can now be entered by pressing the MEMORY button.
9. Press the TIME button again. The last two characters (Minutes) now blink.
10. The exact time can now be entered by pressing the MEMORY button.
11. Press TIME button (or TIME / DATE or TIME): the unit of measurement will flash.
12. Press the "MEMORY to set the unit of measurement (mmHg or kPa)
13. Now after all settings have been made, press the TIME button once again. The date is briefly displayed and then the time. The input is now confirmed, and the clock begins to work.
14. Rapid user select button press the user button, then  is blinking, press again to change user.

To be Noted

With each press of the button (TIME, MEMORY) one input is made (e.g. switching over from hours to minutes mode or altering the value by +1). However, if you keep the respective button pressed, you can switch more quickly to the desired value respectively.

5. CARRYING OUT A MEASUREMENT

5.1. Before the Measurement

To ensure a reliable reading follow these recommendations:

- Avoid eating, drinking alcohol, bathing, smoking as well as any form of exertion directly before taking the measurement as all these factors influence the measurement result.
- Try to rest for 15-30 minutes by sitting in an armchair in a quite atmosphere prior to taking a measurement.
- Always take measurements in a quiet place.
- Stress raises blood pressure. So, don't take measurements during stressful times.
- Measure always on the same arm (preferably left).
- The cuff (arm) should be at the same level as your heart.
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.
- Avoid wearing tight clothing on your arm.
- Do not talk, just remain still during the measurement.

To be Noted

Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.




5.2. Fitting the cuff

Insert air connector into air outlet shown in the given picture and make sure the air connector is fitted properly to avoid air leakage.




Positioning the Cuff on The Arm Correctly


a) The distance between the edge of cuff and the elbow should be approx 2-3cm.



b) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight. Ensure that no space should remain between the cuff and the arm.



c) Lay the arm on a table, with the palm upwards. Support the arm a little with a rest (cushion), so that the cuff rests at about the same height as the heart. Take care, that the cuff lies free. Sit quietly in the same position for 2 minutes after beginning with the measurement.




d) Let legs be uncrossed, feet at on the floor, back and arm supported.



5.3. Start Taking BP Measurement

a) Press the ON/OFF button, the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.





b) After reaching the inflation pressure, the pump stops and the pressure slowly falls away. The cuff-pressure is displ during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and for every pulse beat.

c) When the measurement has been concluded, the measured systolic and diastolic blood- pressure values as well as the pulse frequency are shown on its display.


d) The measurement results are displayed, until you switch the device off. If no button is pressed for 3 minutes, the device switches automatically off, to save the batteries.



Example (Fig.): Systole 118,
Diastole 73, Pulse 75


5.4. Discontinuing a Measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the "ON/OFF" power button can be pressed at any time. The device then immediately lowers the cuff-pressure automatically.



5.5. Memory-Storage and Recall of The Measurements

The blood-pressure monitor automatically stores each of the last 120 measurement values with user identity. By pressing the MEMORY button, an average value of the last 3 measurements as well as the last measurement and the further last 120 measurements (MR119, MR118, ..., MR1) can be displayed one after the other.




(MR1: Values of the last measurement) (MR2-MR120: Values of the measurement before MR1)



5.6. Memory full

Pay attention that the maximum memory capacity is not exceeded. When the memory is full, the old values are automatically deleted. Overwritten with new ones. When memory is full, the display shows the following picture for 1 sec to remind you memory is full.



5.7. Memory– Cancellation of All Measurements

Before you delete all readings stored in the memory, make sure you don't need refer them for future reference. Keeping a written record is prudent and may provide additional information for your doctor's visit. In order to delete all stored readings, press the MEMORY button for at least 5 seconds, the display will show the symbol «CL» and then release the button. To clear the memory permanently, press the MEMORY button while «CL» is flashing to indicate deletion of stored data.



6. Appearance of the Heart Arrhythmia Indicator for Early Detection

This symbol indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal blood pressure-repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times in a week with measurements taken daily) we advise you to tell your doctor.

7. ERROR MESSAGES /MALFUNCTIONS

If an error occurs during a measurement, the measurement is discontinued, and a corresponding error code is displayed.

Error No.	Possible cause(s)
ERR 1	No pulse has been detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the Measurement (Artifact).
ERR 3	The inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 4	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take other reading following following directions carefully. Contact you doctor if you continue to get unusual readings.
ERR 5	Pressure in cuff is over 290mmHg



Other Possible Malfunctions And their Elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains empty when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none">1. Check batteries for correct polarity and if necessary, insert correctly.2. If the display is unusual, re-insert batteries or exchange them
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none">1. Check the positioning of the cuff.2. Measure the blood-pressure again in peace and quiet under observance of the details made under point 5.

Every measurement produces a different value although the instrument functions normally and the values displayed are normal.

1. Please read the following information and the points listed under «Common sources of error». Repeat the measurement.
2. **Please Note:** Blood pressure fluctuates continually so successive measurements will show some variability

Blood pressure measured differs from those values measured by the doctor.

1. Record the daily development of the values and consult your doctor.
2. **Please Note:** Individuals visiting their doctor frequently experience anxiety which can result in a higher reading at the doctor than obtained at home under resting conditions

To Be Noted For licensing, the device has been subjected to strict clinical tests, by which the computer program used to measure the blood-pressure values was tested by experienced specialist doctors in Germany. The same computer program is used in every individual device and has thus also been clinically tested. The manufacture of the devices takes place according to the terms of the European standard for blood-pressure measuring devices (see technical data) you must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. Never attempt to repair the instrument yourself!



8. CARE AND MAINTENANCE, RE-CALIBRATION

- a) Do not expose the device to extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with damp cloth and soapsuds. The cuff must not be washed!
- d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

Periodical Re -calibration Sensitive measuring devices must be checked from time to time for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years.

9. BATTERY LIFE

1000 times measurement with 4- size " AA " alkaline Batteries

10. SAFETY, CARE AND DISPOSAL

Safety & Protection

- This instrument may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for the damage caused by incorrect application.
- This instrument comprises of sensitive components and must be treated with caution. Observe the storage and operating condition described in the "Technical specifications" section!
- Protect it from
 - water and moisture
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - heat and cold
 - direct sunlight
- The cuffs are sensitive and must be handled with care.
- Only pump up the cuff once fitted.
- Do not use the instrument close to strong electromagnetic fields such as mobile telephones or radio installations.
- Do not use the instrument if you think it is damaged or notice anything unusual.
- If the instrument is not going to be used for a prolonged period, the batteries should be removed.
- Read the additional safety instructions in the individual sections of this booklet.
- Ensure that children do not use the instrument unsupervised as some parts are small enough to be swallowed.
- A warning to remove primary batteries if the instruments is not likely to be used for some time.
- Must use the recognized accessories, detachable parts and materials, if the use of other parts or materials can degrade minimum safety.



Warning

1. No modification of this equipment is allowed.
2. Do not modify this equipment without authorization of the manufacturer
3. If this equipment is modified, appropriate inspection must be conducted to ensure continued safe use of the equipment.

Instrument care

Clean the instrument only with a soft, dry cloth

Disposal



Batteries and electronic instruments must be disposed off in accordance with the locally applicable regulations, not with domestic waste.

11. REFERENCE TO STANDARDS

Device standard: Device corresponds to the requirements of the European standard for non-invasive blood pressure monitor

IEC60601-1-6:2010+A1:2013/ EN60601-1-6:2010+A1:2015

IEC60601-1:2005+A1:2012/EN60601-1:2006+A11:2011+A1:2013+A12:2014

IEC60601-1-2:2014/ EN60601-1-2:2015

IEC/EN60601-1-11:2015

IEC80601-2-30:2009+A1:2013/EN80601-2-30:2010+A1:2015

The stipulations of the EU-Guidelines 93/42/EEC for Medical Products Class IIa have been fulfilled.

12. TECHNICAL SPECIFICATIONS







Measurement Procedure:	Oscillometric, corresponding to Korotkoff method: Phase I: systolic, Phase V: diastolic
Display:	Digital display
Measuring range:	Pressure: 30 to 280 mmHg (in 1 mmHg increment) Pulse: 40 to 199 beat/minute
Static accuracy:	Pressure: ± 3 mmHg / Pulse: $\pm 5\%$ of reading
Measuring resolution:	1mmHg
Inflation:	Automatic inflation by internal pump
Memory function:	2 x 120 memories for 2 users (SYS, DIA, Pulse)
Decompression:	Constant exhaust valve system
Power source:	4- size " AA " alkaline Batteries
Rated voltage:	Micro USB DC5.0V/1.0A (direct current)
Operation temperature:	5~40°C/41~104°F
Operation humidity:	15%~85%RH maximum
Storage temperature:	-10~+55°C/14~+131°F
Storage humidity:	10%~95%RH maximum
Dimensions:	142 x 98 x 53 mm



Weight:	494 g (including batteries and cuff)
Cuff pressure display range:	0~290mmHg/0~38.7KPa
Electrical shock protection:	Internal power unit
Safety classifications:	Type BF equipment
Mode of operation:	Continuous operation
Protection against ingress of water:	IP22
Accessories:	M-L size Cuff, 4 "AA" batteries, instruction manual

Please be noticed the power adapter is not supplied from the origin, users can buy the adapter in the market which must comply to EN60601-1,EN60601-1-2

13. REMARK

	Some electrical and electrical equipment forbid to abandon and disposal at will		Reading Instruction Book before use
	Manufacturer's name and address		Type BF equipment
	Inapplicable baby		Attention consult accompanying documents



Cuff Connector



Keep Dry

14. MANUFACTURER'S DECLARATION

The Dr Trust BP Smart is intended for use in the electromagnetic environment specified below. The customer or the user of the Dr Trust BP Smart should assure that it is used in such an environment.

Electromagnetic Emissions: (IEC60601-1-2)

Emission Test	Compliance	Electromagnetic Environment
RF emission CISPR 11	Group 1	The Dr Trust BP Smart uses RF energy only for internal functions. Therefore, this RF emission is extremely weak and there is little chance of it creating any kind of interference whatsoever with nearby electronic equipment.




Emission Test	Compliance	Electromagnetic Environment	
RF emissions CISPR 11	Class B	The Dr Trust BP Smart is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Not applicable		
Voltage fluctuations/flicker IEC 61000-3-3	Not applicable		
Electromagnetic Immunity: (IEC60601-1-2)			
Immunity test	IEC60601-1-2 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electric fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (95% dip in UT.) for 0.5 cycle 40 % UT (60% dip in UT)for 5 cycles 70 % UT (30% dip in UT) for 25 cycles <5 % UT (95% dip in UT) for 5 sec.	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the upper arm style requires continued operation during power mains interruptions, it is recommended that the Dr Trust BP Smart be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Not applicable	Not applicable

Note: UT is the a.c. mains voltage prior to application of the test level.



Immunity test	IEC60601-1-2 test level	IEC60601-1-2 test level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 80% AM (2Hz)	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Dr Trust BP Smart including cables, than their commended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommend separation distance 3V $d = 1.2 \times p^{1/2}$ 80MHz to 800 MHz $d = 2.3 \times p^{1/2}$ MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

Radiated RF IEC 61000-4-3	3 Vrms 80 MHz to 2.5 GHz 80% AM (2Hz)	3 V/m	and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
<p>Note1: At 80 MHz and 800 MHz, the higher frequency range applies. Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			



- a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Dr Trust BP Smart is used exceeds the applicable RF compliance level above, the Dr Trust BP Smart should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Dr Trust BP Smart.
- b Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

Recommended Separation Distances:

Recommended separation distance between portable and mobile RF communications equipment and the Dr Trust BP Smart

The Dr Trust BP Smart is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Dr Trust BP Smart can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Dr Trust BP Smart as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2 \times p^{1/2}$	80 MHz to 800 MHz $d = 1.2 \times p^{1/2}$	800 MHz to 2.5 GHz $d = 2.3 \times p^{1/2}$

0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



● CUSTOMER SUPPORT ●

CONTACT ADDRESS

USA

Nureca INC.USA

276 5th Avenue, Suite 704-397,

New York (NY)-10001, USA

INDIA

Corporate Office (Mumbai)

Nureca Limited

Office No. 101, 1st Floor, Udyog Bhavan,

Sonawala Lane, Goregaon East,

Mumbai, Maharashtra - 400063.

Mfg. By :

Nureca Technologies Pvt. Ltd.

Plot No. 99, 100, Sector 82, JLPL, Mohali, SAS Nagar, Punjab - 160055

Tel. : +91 0172-5292930, 7527013265, 9356658436

Contact us

India : +91-7527013265 /+91-9356658436

Website : www.drtrustusa.com

Crop Website : www.nureca.com

Email : customercare@nureca.com

Connect with us on social networks

Facebook : @drtrust Instagram : @drtrustisin Youtube : NurecaUsa

COPYRIGHT©2021 DR TRUST ALL RIGHTS RESERVED

AR009-00-01



Scan to View
Product Demo Video
www.drtrustusa.com/101