

Control D PRO NEBULIZER



Specifications

Power: 230 V AC, 50 Hz

Maximum Power : 75 W

Medication Capacity : 5 mL

Particle Size: 0.4 to 8 μm

MMAD: 4 μm

Sound Level: 70 dBA(1m)

Weight: 1.4 KG

Nebulization Rate: 0.2 mL or above

Dimensions: L 190 x W 175 x H 105 mm

Compressor Pressure Range: 28 PSI

Operating Pressure Range: 8 to 16 PSI

Compressor Air Output: 10 Lpm

Operating Temperature Range: 0 C to 40 C

Operating Humidity Range: 10 to 75% RH

Storage Humidity Range: 90% RH (Max)

Control D PRO NEBULIZER

PRO NEBULIZER WITH HANDLE & Adjustable Medicine Chamber

Mask Kit Mask Kit includes Adult Mask, Child Mask, Mouth Piece, Air Tube, Medicine Chamber, Air Filters

