

# Control D RESPOCARE NEBULIZER



## Specifications

Power: 230 V AC, 50 Hz  
Maximum Power: 75 W  
Medication Capacity: 5 mL  
Particle Size: 0.4 to 8  $\mu\text{m}$   
MMAD: 4  $\mu\text{m}$   
Sound Level: 70 dBA (1m)  
Weight: 1.4 KG

Nebulization Rate 0.2 mL or above  
Dimensions L 190 x W 175 x H 105 mm  
Compressor Pressure Range 28 PSI  
Operating Pressure Range 8 to 16 PSI  
Compressor Air Output 10 Lpm  
Operating Temperature Range 0°C to 40°C  
Operating Humidity Range 10 to 75% RH  
Storage Humidity Range 90% RH (Max)

## Contents

Nebulizer 1N, Child Mask 1N, Adult Mask 1N, Medicine Chamber 1N, Air Tube 1N, Air Filters 4N

HAIDEN TECHNOLOGY PRIVATE LIMITED, Khasra no 13/4 13/5, Plot No. 8, Gali no. 14, Railway Station Road, Samaipur, New Delhi (North West) – 110042, INDIA