



# NEBULIZER HANDYNEB **GOLD**



**PATENTED**  
DESIGN **MADE IN INDIA**

## KEY FEATURES



Piston Type Compressor



Extensive Ventilation



Safety Fuse - Better Safety



2 Part Medicine Bottle (Integrated Baffle)



Low Noise Level - Less Than 60dBA



Wall Mountable Design



Anti-Shock Case



3 Ball Bearings For Smoother Running

## SAHI CHUNO, **MEDTECH** LO.

### MEDTECH LIFE PVT. LTD.

Regd. Address: B6, Byculla Service Industries, D. K. Marg, Byculla, Mumbai - 400027, Maharashtra, India.

For India	+91 82917 26884	sales@medtechlife.com
For Export	+91 91369 23382	info@medtechlife.com

