



BM 35

beurer  
medical

## BM 35 – Blood pressure monitor

 german  
engineering

risk  
indicator



Risk indicator:  
Classification of measurements with coloured scale

irregular  
heart rhythm



Arrhythmia detection: Warning in case of possible  
heart rhythm disturbance



Fully automatic blood pressure and pulse measurement on the upper arm

Two user memories for 60 measured values each

With cuff for upper arm circumferences from 22-36 cm

Average value of all stored measurements and morning and evening blood  
pressure for the last 7 days

Clear, easy-to-read display

Date and time

Automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Incl. 4 x 1,5V AAA batteries and storage pouch

Medical device

Product measurements (LxWxH): 135 x 105 x 53 mm

Product weight: approx. 326 g (without batteries, with cuff)

XL cuff (30-42 cm) available separately: Art. No. 162.973

5 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

EAN no.: 4211125 65409 2

Item no.: 654.09