

3M™ LITTMANN®

Cardiology IV™

STETHOSCOPE

The sounds you
barely hear are
often the most
important.



Explore the Cardiology IV stethoscope
at www.littmann.co.uk

3M

Littmann®
Stethoscopes

3M™ LITTMANN®

Cardiology IV™

STETHOSCOPE

Hear a new level of quiet.

Catch more of the subtle sounds that are critical to good assessments and diagnoses with the Cardiology IV stethoscope. Improved acoustics and a new design make it more versatile than ever.

Buy now.

3M

3M United Kingdom PLC
Charnwood Campus
10 Bakewell Road
Loughborough, Leicestershire, LE11 5RB
www.littmann.co.uk



Better audibility of high-frequency sounds*



Tunable diaphragms that let you hear both high and low frequency sounds without turning over the chestpiece



Adult and paediatric sides for more versatility



Single-piece diaphragm that's easy to attach and clean because the surface is smooth without crevices

** Adult side, firm pressure. Compared to the 3M™ Littmann® Cardiology III™ Stethoscope.*

3M, Littmann, the Littmann Quality seal, the L, Cardiology III and Cardiology IV are trademarks of 3M. Please recycle. © 3M 2018. All rights reserved.