

## minimus® II

Exact precision that lasts for a lifetime.



 **BIHS**



Rugged, easily operated, in pioneering quality. The chromium-plated screwed tube connection at the top and the ergonomic spoon-shaped grip ensure that **minimus® II** can be operated without hindrance at any time. Also optionally available in a latex-free version (when ordering enter LF for latex-free).

- › Compact, single-tube, graduated-scale sphygmomanometer
- › Highly polished chromium-plated, rugged metal housing with metal ring, corrosion-protected
- › Latex bulb with chromium-plated air release valve
- › Metal valve base
- › Precision air release valve with fine adjustment and wear-free
- › Microfilter protects air release valve and measuring system
- › Spoon-shaped grip of stainless steel
- › Metal connector for simple and fast changes of cuffs
- › Maximum error tolerance  $\pm 3$  mmHg.
- › Linear aluminum scale  $\varnothing$  49 mm, easily legible up to 300 mmHg
- › In vinyl zipper case
- › Comprehensive selection of different sizes and versions of cuffs
- › Specially tempered copper-beryllium diaphragm almost non-aging
- › Approved by the BIHS validation

<b>1312-100</b>	minimus® II, manometer, w/o cuff	<b>1312</b>	minimus® II, velcro cuff (adult)
<b>1310</b>	minimus® II, bandage cuff (adult)	<b>1312-122</b>	minimus® II, velcro cuff (large adult)
<b>1311</b>	minimus® II, hook cuff (adult)	<b>1312-129</b>	minimus® II, velcro cuff (children)

Please visit [riester.de](http://riester.de) for more information