

## minimus® III

Exact precision that lasts for a lifetime.



 **BIHS**

The embodiment of palm-style blood-pressure measuring technology!

Designed for precision and longevity, **minimus® III** is provided with a twin tube, a chromium-plated metal housing and a metal valve base. Also optionally available in a latex-free version (when ordering enter LF for latex-free).

- › Latex bulb with chromium-plated air release valve
- › Precision air release valve with fine adjustment and wear-free
- › Microfilter protects air release valve and measuring system
- › Spoon-shaped grip of stainless steel
- › Specially tempered copper-beryllium diaphragm almost non aging
- › Metal connector for simple and fast changes of cuffs
- › Maximum error tolerance  $\pm 3$  mmHg.
- › Linear aluminum scale  $\varnothing$  49 mm, easily legible up to 300 mmHg
- › In vinyl zipper case
- › Comprehensive selection of different sizes and versions of cuffs
- › Approved by the BIHS validation

<b>1342-100</b>	minimus® III, manometer, w/o cuff	<b>1342</b>	minimus® III, velcro cuff (adult)
<b>1340</b>	minimus® III, bandage cuff (adult)	<b>1342-123</b>	minimus® III, velcro cuff (large adult)
<b>1341</b>	minimus® III, hook cuff (adult)	<b>1342-130</b>	minimus® III, velcro cuff (children)

Please visit [riester.de](http://riester.de) for more information