



beurer
health and well-being



Our products 2024

INTERNATIONAL

www.globalbiomedicals.com

+91-982-8122-223

Health

Wellbeing

Home

Personal Care

Beauty

Fitness

Baby



German
Engineering

Blood pressure and ECG



Detect signs of a stroke before it's too late.

The **BM 53** and **BM 64** Upper Arm Blood Pressure Monitors have an integrated **AFib** function that detects **atrial fibrillation** to help prevent strokes.



Blood pressure



Upper Arm Blood Pressure Monitor BM 53



Upper Arm Blood Pressure Monitor BM 64



Triple measurement

Recommended by doctors for even better accuracy

Made in Europe

Beurer quality standard – made in Europe



Clinically validated

Reliable measurement, even during pregnancy



Bluetooth

App connection

Keep an eye on your blood pressure values at all times with the app "beurer HealthManager Pro"



Inflation technology

Ultimate measurement comfort thanks to gentle build-up of pressure and fast measurement



Advanced arrhythmia detection –

detects atrial fibrillation¹⁾ to help prevent strokes

Atrial fibrillation is the most common form of heart arrhythmia. People affected by it often have no symptoms, but they are at higher risk of a stroke.

A blood pressure monitor with an integrated function to **detect atrial fibrillation (AFib)** is the first step in controlling and monitoring your blood pressure at home. When measuring your blood pressure, the **blood pressure monitor checks how regular your pulse is**. If AFib (atrial fibrillation) is identified, you should always consult a doctor for further advice.

To obtain **even more precise measurements** at home, use the **triple measurement** function. Cardiologists and hypertension societies recommend measuring your blood pressure several times in a row, even when measuring at home.

Don't take any chances when it comes to your heart. We're here to help.

¹⁾ Also known as AFib.

+91-982-8122-223



Blood Pressure Monitor BM 59 – hassle-free measurement without tubes or cables.

Rechargeable battery

Deflation technology

The measurement is taken during deflation

Simple and precise measurement

No tubes, no cables




App connection

With the app "beurer HealthManager Pro"¹⁾



"beurer MyHeart" – a 30-day programme to perfectly complement your blood pressure monitor.

The holistic "beurer MyHeart" concept in the app "beurer HealthManager Pro" helps you integrate a healthier lifestyle into your everyday routine. The four elements – **recipes, exercise, daily inspiration and useful information** – piece together to help you live more healthily. A "beurer MyHeart" health coach can also be booked via the app as an additional consultation service for a fee.




"beurer MyHeart"

Digital coach for advice and guidance

4 building blocks for a healthy daily routine

Nutrition, exercise, useful information and inspiration



Inspiration for greater mindfulness

Motivation that keeps going

Delicious recipes for your well-being

Easy to follow, simple to make and delicious



The new app feature for your digital blood pressure diary.



"Scan & Save" – keep track of your daily measurements in the app "beurer HealthManager Pro". Quick and easy.

Writing down your blood pressure values is a thing of the past! Thanks to the "Scan & Save" function, you can now **save measurements from all Beurer blood pressure monitors digitally in the app** – even without *Bluetooth*[®]. All you need to do is **scan the measurements with a smartphone camera**. The result is automatically detected and saved in the app "beurer HealthManager Pro".

It's easy:

1

Download the App
"beurer HealthManager Pro"



2

Open the "Scan & Save" function in the app and scan the measurements



3

The measurements are automatically saved in the app "beurer HealthManager Pro"



Discover the world of Beurer apps.

We provide devices as well as matching apps to assist you with your day-to-day life.

Scan the QR code to find out whether your smartphone, tablet or laptop is compatible with the Beurer apps. It will show you the system requirements and compatible devices for each app and software program.





Should you measure your blood pressure on your upper arm or wrist? It's up to you.

You can measure your blood pressure **either on your upper arm or on your wrist**. As various tests show, the measurements are **just as reliable either way**.

Upper arm blood pressure monitors

- The **upper arm cuff** helps ensure the **correct position at heart level for an accurate reading**
- An upper arm device is recommended for the elderly, smokers and patients with cardiac arrhythmia and diabetes
- Take a look at our **video tip** for correct measurement on the upper arm:



Wrist blood pressure monitors

- **Simple and easy to use** – handy, compact and ideal when traveling
- To obtain a reliable result, it is important that the wrist is at heart level during the entire measurement
- Take a look at our **video tip** for correct measurement on the wrist:



Upper arm blood pressure monitors

Classic Line

			
BM 26 Item no.: 652.28	BM 28 Item no.: 658.13	BM 35 Item no.: 654.09	BM 40 Item no.: 658.15

Comfort Line

			
BM 44 Item no.: 655.01	BM 45 Item no.: 658.06	BM 49 Item no.: 656.28, item no.: 656.31 Item no.: 656.32, item no.: 656.40	BM 51 easyClip Item no.: 652.04
		 NEW 01/2024	
BM 54 Item no.: 655.12	BM 58 Item no.: 655.16	BM 59 Item no.: 102.27	

Premium Line

	 NEW	 NEW 01/2024		
BM 53 Item no.: 654.59	BM 64 Item no.: 654.00	BM 81 easyLock Item no.: 655.11	BM 85 Item no.: 658.03	BM 96 Cardio Item no.: 658.01

Wrist blood pressure monitors

Classic Line

		
BC 27 Item no.: 659.04	BC 28 Item no.: 650.44	BC 30 Item no.: 650.51

Comfort Line

		
BC 21 Item no.: 650.46	BC 44 Item no.: 659.05	BC 51 Item no.: 650.60

Premium Line

	
BC 54 Item no.: 650.54	BC 87 Item no.: 657.03

BM 26



Upper arm blood pressure monitor

- Risk indicator: Classification of measurements with colour scale
- Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the upper arm

Average of the last 3 measurements

Clear, easy-to-read display

Four user memories for 30 measurements each

Risk indicator

Arrhythmia detection

Medical device

For upper arm circumferences from 22 – 35 cm

Date and time/automatic switch-off

Incorrect usage message

Low battery indicator

Storage bag

Incl. 4 x 1.5V AA batteries

XL cuff (30 – 42 cm) Item no.: 162.973 available as an accessory

5-year guarantee

Pack: 4 Shipping carton: 6 x 4

EAN no.: 4211125652289

Item no.: 652.28



BM 28



Upper arm blood pressure monitor

- Resting indicator¹⁾ for accurate measurements
- Cuff position control: For ensuring the cuff is correctly positioned on the upper arm



Fully automatic blood pressure and pulse measurement on the upper arm

Cuff position control

Measurement reminder: Set 2 alarm times

Average of all saved measurements

Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read display

Four user memories for 30 measurements each

Risk indicator

Arrhythmia detection

Medical device

For upper arm circumferences from 22 – 42 cm

Date and time/automatic switch-off

Incorrect usage message

Low battery indicator

Storage bag

Incl. 4 x 1.5V AA batteries

Mains adapter available separately Item no.: 071.95

BM 28 On-Pack (= incl. mains adapter)

EAN no.: 4211125658304 Item no.: 658.30

5-year guarantee

Pack: 4 Shipping carton: 4 x 4

EAN no.: 4211125658137

Item no.: 658.13

¹⁾ Patent EP 1 673 009

BM 28 product info

Reliable measurement, even during pregnancy

It is advised to regularly monitor your blood pressure during pregnancy – especially for those with high blood pressure. The BM 28 has been clinically validated for use in pregnancy, making it your ideal companion for optimal pregnancy care.



Before, during or after pregnancy, women may experience high blood pressure, which increases the risk of pre-eclampsia.

For this reason, **regular and, above all, reliable blood pressure measurements are extremely important for you and your baby.**

In **pregnant women with high blood pressure**, changes in the tissue and artery structures triggered by pregnancy can mean that **standard blood pressure monitors display deviating measured values.**

A clinical validation study of the **BM 28 Oscillometric Blood Pressure Monitor from Beurer confirms the measurement accuracy** of the blood pressure values in pregnant women and also in cases of existing arterial hypertension and pre-eclampsia.



Extra tip: With its overall technology, the BM 28 Upper Arm Blood Pressure Monitor is also suitable for regular monitoring of blood pressure at home for those who are not pregnant.



BM 35



Upper arm blood pressure monitor

- Risk indicator: Classification of measurements with colour scale
- Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the upper arm

Average of all saved measurements

Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read display

Two user memories for 60 measurements each

Risk indicator

Arrhythmia detection

Medical device

For upper arm circumferences from 22 – 36 cm

Date and time/automatic switch-off

Incorrect usage message

Low battery indicator

Storage bag

Incl. 4 x 1.5V AAA batteries

XL cuff (30 – 42 cm) Item no.: 162.973 available separately as an accessory

5-year guarantee

Pack: 4 Shipping carton: 4 x 4

EAN no.: 4211125654092

Item no.: 654.09

BM 40



Upper arm blood pressure monitor

- Risk indicator: Classification of measurements with colour scale
- Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the upper arm

Average of all saved measurements

Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read XL display

Two user memories for 60 measurements each

Risk indicator

Arrhythmia detection

Medical device

For upper arm circumferences from 22 – 35 cm

Date and time/automatic switch-off

Incorrect usage message

Low battery indicator

Storage bag

Incl. 4 x 1.5V AA batteries

XL cuff (30 – 42 cm) Item no.: 162.973 available separately as an accessory

Mains adapter available separately Item no.: 071.95

BM 40 On-Pack (= incl. mains adapter)

EAN no.: 4211125658175 Item no.: 658.17

5-year guarantee

Pack: 4 Shipping carton: 3 x 4

EAN no.: 4211125658151

Item no.: 658.15

BM 44



Upper arm blood pressure monitor

- Easy to use: Single-button operation
- Blue illuminated, easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the upper arm

Illuminated START/STOP button

Blue illuminated, easy-to-read XL display

Risk indicator

Arrhythmia detection

Medical device

For upper arm circumferences from 22 – 30 cm

Automatic switch-off

Incorrect usage message

Low battery indicator

Storage bag

Incl. 4 x 1.5V AA batteries

XL cuff (30 – 42 cm) Item no. 162.795 available separately as an accessory

Mains adapter available separately Item no. 071.95

5-year guarantee

Pack: 4 Shipping carton: 3 x 4

EAN no.: 4211125655013

Item no.: 655.01

BM 45



Upper arm blood pressure monitor

- Risk indicator: Classification of measurements with colour scale
- White illuminated, easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the upper arm
Easy to use via touch sensor buttons

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Illuminated START/STOP button
White illuminated, easy-to-read XL display
Two user memories for 60 measurements each
Risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 36 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AA batteries
XL cuff (35 – 44 cm) Item no. 163.387 available separately as an accessory
Mains adapter available separately Item no. 071.95

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125658069
Item no.: 658.06

BM 49



Speaking upper arm blood pressure monitor

- With voice output of measurements in multiple languages
- Blue illuminated, easy-to-read XL display
- Adjustable volume



Fully automatic blood pressure and pulse measurement on the upper arm
Voice output:
(a) DE, FR, IT, NL (b) EN, ES, RU, GR
(c) RO, PL, CS, HU (d) FA, AR, EN

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Blue illuminated, easy-to-read XL display
Two user memories for 60 measurements each
Risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 36 cm
Date and time/automatic switch-off
Incorrect usage message
Storage bag
Incl. 4 x 1.5V AA batteries
XL cuff (35 – 44 cm) Item no. 163.387 available separately as an accessory
Mains adapter available separately Item no. 071.95
5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: (a) 4211125656287
Item no.: (a) 656.28
EAN no.: (b) 4211125656409
Item no.: (b) 656.40
EAN no.: (c) 4211125656317
Item no.: (c) 656.31
EAN no.: (d) 4211125656324
Item no.: (d) 656.32

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125656324
Item no.: (d) 656.32

BM 51 easyClip



Upper arm blood pressure monitor

- Innovative easyClip cuff (22 – 42 cm): Cuff is easy to fit with just one action
- Inflation technology: Fast and convenient measurement that starts while the cuff is inflating
- White illuminated, easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the upper arm
Automatic measurement during inflation: For quick and convenient measuring

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
White illuminated XL display
Two user memories for 100 measurements each
Risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 42 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AA batteries
Mains adapter available separately Item no.: 071.95

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125652043
Item no.: 652.04



BM 53



Upper arm blood pressure monitor

- Advanced arrhythmia detection: Detects atrial fibrillation to help prevent strokes
- Triple measurement is recommended by doctors for even more accurate measurements
- Reliable measurement, even during pregnancy



Fully automatic blood pressure and pulse measurement on the upper arm
Resting indicator¹⁾ for accurate measurements
Cuff position control

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Automatic pressure preselection and deflation
Very easy-to-read XL display with backlight
Two user memories for 100 measurements each
Risk indicator
Arrhythmia detection
Medical device
Premium cuff for upper arm circumferences from 22 – 42 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AAA batteries
Mains adapter available separately

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125654597
Item no.: 654.59

BM 54



Upper arm blood pressure monitor

- Universal cuff, including for larger upper arm circumferences up to 44 cm
- Risk indicator: Classification of measurements with colour scale



Fully automatic blood pressure and pulse measurement on the upper arm

Wireless transfer of measurements
Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Clear, easy-to-read XL display
Two user memories for 60 measurements each
Risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 44 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AAA batteries
Free app "beurer HealthManager Pro"

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125655129
Item no.: 655.12

BM 58



Upper arm blood pressure monitor

- Transfer measurements to the app via USB using the "beurer HealthManager Pro" USB uploader
- Easy to use via touch sensor buttons
- Clear, easy-to-read, black XL display



Fully automatic blood pressure and pulse measurement on the upper arm

PC interface / USB cable included
Cuff holder
Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Clear, easy-to-read, black XL display
Two user memories for 60 measurements each
Risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 30 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AA batteries
XL cuff (30 – 42 cm) Item no. 163.246 available as an accessory
Mains adapter available separately Item no. 071.95
Free app "beurer HealthManager Pro" and USB uploader


5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125655167
Item no.: 655.16

BM 59



Upper arm blood pressure monitor

- Upper arm blood pressure monitor – cordless & automatic measurement on the upper arm
- Deflation technology
- No tubes, no cables – just simple and precise measurement
- Rechargeable battery



Fully automatic blood pressure and pulse measurement on the upper arm
No tubes, no cables – just simple and precise measurement

Wireless transfer of measurements
Rechargeable battery
USB-C port

LED risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 42 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Incl. 4 x 1.5V AAA batteries
Incl. USB-C cable for charging the battery
Free app "beurer HealthManager Pro"


5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125102272
Item no.: 102.27

BM 64



Upper arm blood pressure monitor

- Advanced arrhythmia detection: Detects atrial fibrillation to help prevent strokes
- Triple measurement is recommended by doctors for even more accurate measurements
- Inflation technology: Ultimate measurement comfort thanks to gentle build-up of pressure and fast measurement



Fully automatic blood pressure and pulse measurement on the upper arm
Resting indicator¹⁾ for accurate measurements
Cuff position control

Wireless transfer of measurements
Easy to use via 4 touch sensor buttons
Automatic pressure preselection and deflation

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Very easy-to-read black XL display with backlight
Two user memories for 120 measurements each
Risk indicator
Arrhythmia detection
Medical device
Premium cuff for upper arm circumferences from 22 – 42 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage bag
Incl. 4 x 1.5V AAA batteries
Mains adapter available separately
Free app "beurer HealthManager Pro"

5-year guarantee
Pack: 4 Shipping carton: 4 x 4
EAN no.: 4211125654009
Item no.: 654.00

BM 81 easyLock



Upper arm blood pressure monitor

- Tubular easyLock cuff (24 – 40 cm) for comfortable and easy fitting on the upper arm
- easyLock locking mechanism for a perfect fit and secure hold
- Innovative design: The device is attached directly to the cuff without any tubes



Fully automatic blood pressure and pulse measurement on the upper arm
Inflation technology: Ultimate measurement comfort thanks to gentle build-up of pressure and fast measurement

Wireless transfer of measurements


Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Very easy-to-read black XL display
Two user memories for 120 measurements each
LED risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 24 – 40 cm
Date and time/automatic switch-off
Incorrect usage message
Low battery indicator
Storage box
Incl. 4 x 1.5V AAA batteries
Free app "beurer HealthManager Pro"

5-year guarantee
Pack: 4 Shipping carton: 2 x 4
EAN no.: 4211125655112
Item no.: 655.11



BM 81 easyLock product info

Turn to this intuitive design with integrated cuff and app.



The award-winning **BM 81 easyLock** uses no annoying tubes or cables. Thanks to the **integrated cuff** with **easyLock system**, the device is easy to put on, fits perfectly and is held securely.

Measuring blood pressure has never been this easy or precise! Position the cuff, tighten the easyLock rotary knob and start the measurement. **It doesn't get any better than this!**



1. Put it on
Simple and comfortable thanks to integrated, tubular cuff



2. Tighten
Perfect fit and secure hold – all without a hook-and-loop fastener



3. Take a measurement
Ultimate measurement comfort thanks to gentle build-up of pressure and fast measurement

BM 85



Upper arm blood pressure monitor

- Resting indicator¹⁾ for accurate measurements
- Universal cuff, including for larger upper arm circumferences up to 42 cm
- With lithium-ion battery: Never change batteries again



Fully automatic blood pressure and pulse measurement on the upper arm
Automatic measurement during inflation:
For quick and convenient measuring

Transfer measurements to the app "beurer HealthManager Pro" via Bluetooth® or USB uploader
Wireless transfer of measurements
PC interface / USB power mains adapter²⁾
Cuff holder

Average of all saved measurements
Morning and evening blood pressure averages for the last 7 days
Illuminated START/STOP button
White illuminated, easy-to-read XL display
Two user memories for 60 measurements each
LED risk indicator
Arrhythmia detection
Medical device
For upper arm circumferences from 22 – 42 cm
Date and time/automatic switch-off
Incorrect usage message
Battery charge level
Storage bag
Li-ion battery

Free app "beurer HealthManager Pro" and USB uploader
5-year guarantee
Pack: 4 Shipping carton: 2 x 4
EAN no.: 4211125658038
Item no.: 658.03



BM 96 Cardio



- Blood pressure monitor with ECG function**
- 1-channel ECG for recording heart rhythm with practical ECG stick
 - Optimum monitoring of blood pressure and heart rhythm
 - Inflation technology: Ultimate measurement comfort thanks to gentle build-up of pressure and fast measurement
 - Detects atrial fibrillation (AFib) to help prevent strokes



Fully automatic blood pressure and pulse measurement on the upper arm



- Inflation technology
- Transfer measurements to the app "beurer HealthManager Pro" via Bluetooth® or USB uploader
- ECG recording takes 30-120 seconds using ECG stick
- Morning and evening blood pressure averages for the last 7 days
- Illuminated START/STOP button
- White illuminated, easy-to-read XL display
- Two user memories for 60 measurements each
- Risk indicator
- Medical device
- For upper arm circumferences from 22 - 42 cm
- Alarm function
- Date and time/automatic switch-off
- Incorrect usage message
- Low battery indicator
- Storage bag
- Incl. 4 x 1.5V AA batteries
- Free app "beurer HealthManager Pro" and USB uploader

5-year guarantee
 Pack: 4 Shipping carton: 2 x 4
 EAN no.: 4211125658014
 Item no.: 658.01


"beurer MyCardio Pro" product info

beurer MyCardio Pro – your ECG professional at home

The "beurer MyCardio Pro" additional feature found in the app "beurer HealthManager Pro" provides you with your personal ECG evaluation in just a few seconds. In addition to this analysis, "beurer MyCardio Pro" can also prepare a **detailed medical report**. This report can either be sent directly to your doctor as a PDF file or you can print it out at home for discussion at your next doctor's appointment.



The **ECG measurements** are clearly displayed and **professionally analysed in the app using clinically proven learning algorithms**. Various **health parameters** are included, such as **well-being, stress level and energy value**.



ME 90



- Mobile ECG device**
- Records your heart rhythm
 - Optimum monitoring of the heart rhythm
 - Self-monitoring of heart function at home or on the move
 - Transfer measurements to the app "beurer HealthManager Pro" via Bluetooth® or USB uploader



- For persons with:
- High blood pressure, raised cholesterol values, diabetes, family predisposition
 - Excess weight, lack of physical activity, high consumption of tobacco, stress, lack of sleep
 - History of heart attacks and/or strokes

- Wireless transfer of measurements
- Duration of recording: 30 seconds
- Recording taken from: Chest-arm-hand or hand-arm-upper body-arm-hand
- Immediate information about heart rhythm after completion of recording
- 36 memory spaces
- Medical device
- PC interface / USB cable included
- Date and time/automatic switch-off
- Battery level indicator
- Incl. 2 x 3V CR 2032 lithium batteries
- Free app "beurer HealthManager Pro" and USB uploader. Complete range of uses only available with the app.


5-year guarantee
 Pack: 4 Shipping carton: 8 x 4
 EAN no.: 4211125454159
 Item no.: 454.15

ECG devices product info

ECG recordings made easy – keep an eye on your heart rhythm.

Perform your own ECG recording for the first time with ease. Your heart rhythm is recorded, **displayed simply**, and processed for subsequent interpretation by your doctor.

Our ECG devices **evaluate the measurements** and automatically assess whether the **ECG rhythm** is too fast, too slow, irregular, linked to pauses, or normal. To summarise, an easily intelligible symbol is shown on the displays immediately following the measurement.



An all-in-one solution for your health
 The **BM 96 Cardio** combines blood pressure and pulse monitoring, and ECG recording with one device. The combination of ECG and blood pressure measurement provides the ideal way to comprehensively monitor the health of your heart. This makes the **BM 96 Cardio** the perfect companion for all hypertensive patients who want to oversee and monitor their treatment. It's also a great preventative aid.



BC 21



- Talking wrist blood pressure monitor**
- With voice output of measurements in multiple languages
 - Risk indicator: Classification of measurements with colour scale



Fully automatic blood pressure and pulse measurement on the wrist
 Voice output: DE, EN, FR, IT, TR

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 14 – 19.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 8 x 4
 EAN no.: 4211125650469
 Item no.: 650.46

BC 27



- Wrist blood pressure monitor**
- Risk indicator: Classification of measurements with colour scale
 - Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the wrist

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 14 – 19.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 8 x 4
 EAN no.: 4211125659042
 Item no.: 659.04

BC 28



- Wrist blood pressure monitor**
- Risk indicator: Classification of measurements with colour scale
 - Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the wrist

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 14 – 19.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 8 x 4
 EAN no.: 4211125650445
 Item no.: 650.44

BC 30



- Wrist blood pressure monitor**
- Risk indicator: Classification of measurements with colour scale
 - Arrhythmia detection: Alerts in the case of possible cardiac arrhythmia



Fully automatic blood pressure and pulse measurement on the wrist

Average of the last 3 measurements

Clear, easy-to-read display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 13.5 – 19.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 8 x 4
 EAN no.: 4211125650513
 Item no.: 650.51

BC 44



- Wrist blood pressure monitor**
- Easy to use: Single-button operation
 - Blue illuminated, easy-to-read XL display



Fully automatic blood pressure and pulse measurement on the wrist

Illuminated START/STOP button
 Blue illuminated, easy-to-read XL display


Clear, easy-to-read display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 14 – 19.5 cm
 Automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 4 x 4
 EAN no.: 4211125659059
 Item no.: 659.05

BC 51



- Wrist blood pressure monitor**
- Position indicator: Accurate measurements due to correct measurement position at heart level
 - Inflation technology: Fast and convenient measurement that starts while the cuff is inflating



Fully automatic blood pressure and pulse measurement on the wrist
 Automatic measurement during inflation: For quick and convenient measuring

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Large, very easy-to-read XL display
 Two user memories for 120 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 12.5 – 21 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

5-year guarantee
 Pack: 4 Shipping carton: 4 x 4
 EAN no.: 4211125650605
 Item no.: 650.60



BC 54



Wrist blood pressure monitor

- Inflation technology: Fast and convenient measurement that starts while the cuff is inflating
- Risk indicator: Classification of measurements with colour scale
- Clear, easy-to-read black display



Fully automatic blood pressure and pulse measurement on the wrist
 Automatic measurement during inflation: For quick and convenient measuring
 Wireless transfer of measurements

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Very easy-to-read black display
 Two user memories for 60 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 13.5 – 21.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries

Free app "beurer HealthManager Pro"

5-year guarantee
 Pack: 4 Shipping carton: 4 x 4
 EAN no.: 4211125650544
 Item no.: 650.54

BC 87



Wrist blood pressure monitor

- Wrist blood pressure monitor with patented resting indicator¹⁾ for precise measurement results
- Position indicator shows the optimal measurement position at heart level



Fully automatic blood pressure and pulse measurement on the wrist
 Automatic measurement during inflation: For quick and convenient measuring
 Wireless transfer of measurements

Average of all saved measurements
 Morning and evening blood pressure averages for the last 7 days

Clear, easy-to-read, black XL display
 Two user memories for 120 measurements each
 Risk indicator
 Arrhythmia detection
 Medical device
 For wrist circumferences from 13.5 – 21.5 cm
 Date and time/automatic switch-off
 Incorrect usage message
 Low battery indicator
 Storage box
 Incl. 2 x 1.5V AAA batteries


Free app "beurer HealthManager Pro"

5-year guarantee
 Pack: 4 Shipping carton: 4 x 4
 EAN no.: 4211125657031
 Item no.: 657.03

¹⁾ Patent EP 1 673 009

BC 87 product info


Precise measurements thanks to integrated resting pulse and position check.



The Beurer BC 87 Wrist Blood Pressure Monitor provides a **hassle-free, accurate way** to measure your blood pressure; with **2 intelligent monitoring functions** that make it so simple to use.

1. Resting indicator

The patented resting indicator* monitors your resting pulse during the measurement.



2. Correct position

The position indicator with traffic light function shows at a glance whether the device is positioned correctly at heart level. Now you can rely on accurate measurements when checking your own blood pressure.



Body temperature



Hygienic & safe –
fever measurement without skin contact.

The **infrared technology** of our **non-contact clinical thermometers** allows them to be used without coming into contact with the skin, making them an easy and hygienic way to care for those who need it.



Non-Contact Thermometer FT 100

Measurement of body, ambient and surface temperature

Infrared technology

Non-contact measurement without touching the skin or surface directly

LED fever indicator

Traffic light signal for quick assessment



Reliable


The distance sensor monitors the correct distance to the forehead

Convenient to use

Very easy-to-read XL display, 60 memory spaces, emoticon guidance

Versatile

Measurement of body, ambient and surface temperature



Non-contact clinical thermometers



FT 85
Item no.: 795.06



FT 90
Item no.: 795.31



FT 95 Bluetooth®
Item no.: 795.07




FT 100
Item no.: 795.05



NON-CONTACT MEASUREMENT

Beurer's non-contact clinical thermometers have an **integrated infrared sensor**. This measures the infrared radiation emitted by the body and uses it to determine the body temperature.





FT 09/1 White



Thermometer

- auto off
- 1 memory
- water proof
- incl. protective cover


Measurement time approx. 30 sec.
Contact-measurement technology
Fever indicator: Warning tone from 37.8°C
Display in °C
Waterproof
Can be disinfected
Measurement accuracy: ±0.1°C: 35.5°C - 42.0°C ±0.2°C: 32.0°C - 35.4°C ; 42.1°C - 42.9°C
1 memory space
Automatic switch-off
Acoustic signal when measurement is complete
Medical device
Low battery indicator
With protective cap
Incl. 1 x 1.5V LR41 battery
5-year guarantee
Pack: 12 Shipping carton: 12 x 12
EAN no.: 4211125791155
Item no.: 791.15

BY 11 Monkey / Frog / Dog



Instant thermometer

- Extra-flexible measuring tip



- express 10 seconds
- auto off
- fever indication
- 1 memory
- °C / °F display
- incl. protective cover


Measurement time approx. 10 sec.
Contact-measurement technology
Fever indicator: Warning tone from 37.8°C
Displays measurements in °C and °F
Waterproof tip and display
Can be disinfected
Measurement accuracy: ±0.1°C: 35.5°C - 42.0°C (±0.2°F: 95.9°F - 107.6°F) ±0.2°C: 32.0°C - 35.4°C ; 42.1°C - 42.9°C (±0.4°F: 89.6°F - 95.7°F ; 107.8°F - 109.2°F)
1 memory space
Automatic switch-off
Acoustic signal when measurement is complete
Medical device
With protective cap
Incl. 1 x 1.5V LR41 battery
Incl. sales display
5-year guarantee
Pack: 20 Shipping carton: 10 x 20
EAN no.: Monkey: 4211125950040
Item no.: Monkey: 950.04
EAN no.: Frog: 4211125950057
Item no.: Frog: 950.05
EAN no.: Dog: 4211125950064
Item no.: Dog: 950.06

FT 15/1



Instant thermometer

- Extra-flexible measuring tip



- express 10 seconds
- auto off
- fever indication
- 1 memory
- °C / °F display

Measuring time approx. 10 sec.
Contact-measurement technology
Fever indicator: Warning tone from 37.8°C
Displays measurements in °C and °F
Can be disinfected
Measurement accuracy: ±0.1°C: 35.5°C - 42.0°C (±0.2°F: 95.9°F - 107.6°F) ±0.2°C: 32.0°C - 35.4°C ; 42.1°C - 42.9°C (±0.4°F: 89.6°F - 95.7°F ; 107.8°F - 109.2°F)
1 memory space
Automatic switch-off
Acoustic signal when measurement is complete
Medical device
Low battery indicator
With protective cover
Incl. 1 x 1.5V LR41 battery
Incl. sales display
5-year guarantee
Pack: 20 Shipping carton: 10 x 20
EAN no.: 4211125794101
Item no.: 794.10

FT 58



Ear thermometer

- Hygienic and safe measurement
- With interchangeable protective cap
- Body, object and room temperature
- 10 protective caps included



Measurement in seconds
Infrared measurement
Fever indicator: From 38.0°C, warning tone from 37.5°C
Displays measurements in °C and °F
Object temperature measurement
Room temperature display
Face icons
Smiling face: Everything is OK
Sad face: ≥ 38°C (≥ 100.4°F)
Measurement accuracy: Ear: ±0.2°C: 35.5°C - 42.0°C (±0.4°F: 95.9°F - 107.6°F) Remaining measurement range: From ≤ 35.5°C and ≥ 42°C (≤ 95.5°F and ≥ 107.6°F) ±0.3°C (±0.5°F)
Object: ±4% or ±2°C (±4°F) at -20°C - 100°C (-4°F - 212°F)
Room: ±2°C: 0 - 50°C (±4°F: 32°F - 122°F)
Ejection mechanism
10 memory spaces
Automatic switch-off
Acoustic signal when measurement is complete
Medical device
Low battery indicator
With storage box
Incl. 2 x 1.5V AAA batteries
Incl. 10 protective caps
Set of replacement caps available for purchase (20 caps)
5-year guarantee
Pack: 4 Shipping carton: 12 x 4
EAN no.: 4211125795337
Item no.: 795.33

FT 58 Protective Caps



Set of replacement caps

- 20 protective caps



20 single-use protective caps for the FT 58 Ear Thermometer
For hygiene purposes and to ensure a precise result, replace the protective caps after each measurement
Pack: 10 Shipping carton: 10 x 10
EAN no.: 4211125795344
Item no.: 795.34

FT 65



Multi functional thermometer

- 6-in-1 function: Ear, forehead and surface temperature, fever indicator, date and time, 10 memory spaces



Measurement in seconds
Infrared measurement
LED fever indicator from 38°C
Displays measurements in °C and °F
Can be disinfected
Colour indicator
Green: ≤ 37.9°C (≤ 100.2°F)
Red: ≥ 38°C (≥ 100.4°F)
Measurement accuracy: Ear/forehead: ±0.2°C: 35°C - 42°C (±0.4°F: 95°F - 107.6°F) In remaining measurement range: ±0.3°C (±0.5°F)
Object: ±1.5°C (±2.7°F) at <30°C (86°F), ±5% at ≥30°C (86°F)
10 memory spaces with date and time
Automatic switch-off
Large display
Medical device
Low battery indicator
Incl. 2 x 1.5V AAA batteries
5-year guarantee
Pack: 4 Shipping carton: 8 x 4
EAN no.: 4211125795146
Item no.: 795.14

Pulse oximeters




It's good to know your pulse and your arterial oxygen saturation levels.



Your **heart rate** indicates the rate at which **your heart pumps blood through your body**. Your **pulse** is measured using a pulse oximeter and is an important indicator of your health and fitness level.

Pulse oximetry also measures **arterial oxygen saturation** using two light sources of different wavelengths. As a result, the Beurer pulse oximeters are **perfect for at-risk patients and athletes, both young and old**. Our oximeters are **small, lightweight and fit in any bag**.

Our **new PO 13 Kids model** is adapted to suit **children's fingers**, is **simple to use** and has a fun child-friendly design.



www.globalbiomedicals.com

Pulse oximeters



PO 30
Item no.: 454.30



PO 35
Item no.: 454.31



PO 45
Item no.: 454.34



PO 60 Bluetooth®
Item no.: 454.20



PO 80
Item no.: 454.40

Pulse oximeter with app connection

Pulse oximeter with software connection

Pulse Oximeter PO 13 Kids



Child-friendly

Adapted to suit children's fingers with a friendly giraffe design

2-in-1

Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse)

At a glance

Graphic pulse curve display

Easy and painless

The measurement takes just 8 seconds



High-quality colour display

4 display formats and adjustable display brightness



SUPER SIMPLE

The **very simple and completely painless measurement** takes **only a few seconds** for fuss-free results.



Oxygen keeps your body cells alive – learn more about pulse oximetry.



How pulse oximetry works.


Using pulse oximetry, **arterial oxygen saturation can be measured without having to take blood samples.** The technology measures two light sources of different wavelengths passing through the fingertip.

Our body's oxygen requirements.

The human **body needs approx. 400 litres of oxygen (O₂) each day.** A healthy oxygen saturation level is around 94% to 99%. Checking your level on a regular basis makes it possible to spot the **signs of a deficiency** and take action early.

What is the Perfusion Index?


The Perfusion Index (PI) determines the **blood flow value of the finger.** It can range from 0.3% to 20% and varies depending on the patient, measurement location and state of health.



Oxygen supply.

Oxygen keeps the cells in your body functioning. Transporting a large amount of oxygen requires a "carrier" – **haemoglobin.**

The **oxygen attaches to the haemoglobin** in the red blood cells and is supplied to the organs via the bloodstream.




Who might benefit from a pulse oximeter? Here's what we recommend.

A pulse oximeter is recommended for **athletes** who want to keep an eye on their vital signs. The use of a pulse oximeter is also recommended in situations where the oxygen content of the ambient air is low, e.g. when **you're at higher altitudes for long periods** or during air sports activities.

Measurement of oxygen saturation in the blood also plays an important role for the large number of people who suffer from chronic or acute **respiratory diseases.** Beurer's pulse oximeters are ideal for these patients as they can provide information about problems with the body's oxygen supply.

Our new PO 13 Kids model with a giraffe design has been specially developed for children suffering from heart failure, chronic obstructive pulmonary diseases, bronchial asthma or sleep apnoea.



IMPORTANT.

In high-risk patients, self-monitoring using a pulse oximeter is not sufficient. Medical supervision must also be sought.





PO 13 Kids



- Pulse oximeter**
- Determination of arterial oxygen saturation (SpO₂) and heart rate (pulse frequency)
 - Very easy and completely painless measurement within 8 seconds
 - Displays measurements within seconds without testing the child's patience



Particularly suitable for children with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma, sleep apnoea Adapted to suit children's fingers Easy-to-read colour display Colour display with 4 display formats Adjustable display brightness The display is oriented automatically (portrait, landscape) Graphic pulse display
Small and light – for use at home and on the move
Automatic switch-off Low battery indicator Incl. strap, belt pouch and screwdriver
Incl. 2 x 1.5V AAA batteries Medical device
5-year guarantee Pack: 4 Shipping carton: 4 x 4 EAN no.: 4211125600952 Item no.: 600.95

PO 30



- Pulse oximeter**
- Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse)



Particularly suitable for persons with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma Also suitable for sports at high altitudes Colour display with 4 viewing perspectives Adjustable display brightness Automatic orientation of display Graphic pulse display
Small and light – for use at home and on the move
Automatic switch-off Low battery indicator Incl. strap and belt pouch
Incl. 2 x 1.5V AAA batteries Medical device
5-year guarantee Pack: 4 Shipping carton: 8 x 4 EAN no.: 4211125454302 Item no.: 454.30

PO 35



- Pulse oximeter**
- Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse)



Particularly suitable for persons with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma Also suitable for sports at high altitudes Colour display with 4 viewing perspectives Adjustable display brightness Automatic orientation of display Graphic pulse display
Small and light – for use at home and on the move
Automatic switch-off Low battery indicator Incl. strap and belt pouch
Incl. 2 x 1.5V AAA batteries Medical device
5-year guarantee Pack: 4 Shipping carton: 8 x 4 EAN no.: 4211125454319 Item no.: 454.31

PO 45



- Pulse oximeter**
- Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse) and Perfusion Index (PI)
 - Symbol displayed for inaccurate measurements



Particularly suitable for persons with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma Also suitable for sports at high altitudes Colour display with 7 display formats Adjustable display brightness Graphic pulse display Symbol displayed for inaccurate measurements Perfusion Index Small and light – for use at home and on the move
Automatic switch-off Low battery indicator Incl. strap and belt pouch
Incl. 2 x 1.5V AAA batteries Medical device
5-year guarantee Pack: 4 Shipping carton: 8 x 4 EAN no.: 4211125454340 Item no.: 454.34

PO 60 Bluetooth®



- Pulse oximeter**
- Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse)
 - Transfer measurements to the app "beurer HealthManager Pro" via Bluetooth®




Particularly suitable for persons with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma Also suitable for sports at high altitudes Colour display with 4 viewing perspectives Adjustable display brightness Automatic orientation of display Graphic pulse display Wireless transfer of measurements Small and light – for use at home and on the move Memory space for 100 measurements based on the "finger-in-finger-out" principle Automatic switch-off Low battery indicator Incl. strap and belt pouch
Incl. 2 x 1.5V AAA batteries Medical device Complete range of uses only available with the free app "beurer HealthManager Pro"
5-year guarantee Pack: 4 Shipping carton: 8 x 4 EAN no.: 4211125454203 Item no.: 454.20

PO 80



- Pulse oximeter**
- Measurement of arterial oxygen saturation (SpO₂) and heart rate (pulse)
 - Lithium-ion battery – can be charged by USB
 - USB cable for transferring measurements to your PC



Particularly suitable for persons with: Heart failure Chronic obstructive pulmonary diseases Bronchial asthma Also suitable for sports at high altitudes XL colour display with 4 viewing perspectives Adjustable display brightness Automatic orientation of display Graphic pulse display Optional pulse tone Small and light – for use at home and on the move Signal function – limit values can be set individually Recording up to 24 hours PC software for detailed display Automatic switch-off Battery level display Incl. strap, belt pouch, USB cable, USB charger and software Metal storage box Incl. lithium-ion battery: Can be charged by USB Medical device Free "SpO ₂ Assistant" PC software
5-year guarantee Pack: 4 Shipping carton: 8 x 4 EAN no.: 4211125454401 Item no.: 454.40