

# 3M™ Littmann® Classic II Pediatric and Littmann Classic II Infant Stethoscopes



3M Medical  
*from the 3M Health Care family*



# 3M™ Littmann®

## Classic II Pediatric and Classic II Infant Stethoscopes



Two-sided models with traditional bell/diaphragm chestpieces are designed, sized and acoustically precise for children and infants.

### Features

- One-inch traditional bell (pediatric) and ¾" bell (infant) combined with a floating diaphragm
- Solid stainless steel chestpiece
- Excellent acoustic seal and comfortable fit with patented 3M™ Littmann® Snap Tight Soft-Sealing eartips
- Comfortably angled, anatomically correct headset
- Single-lumen tubing available in attractive colors
- Three-year warranty

### 3M Littmann Classic II Pediatric Stethoscope Ordering Information

Catalog #	Length	Tube Color
2113	28 in. (71cm)	Black
2113R	28 in. (71cm)	Red
2119	28 in. (71cm)	Caribbean Blue
2122	28 in. (71cm)	Raspberry
2123	28 in. (71cm)	Navy Blue

### 3M Littmann Classic II Infant Stethoscope Ordering Information

Catalog #	Length	Tube Color
2114	28 in. (71cm)	Black
2114R	28 in. (71cm)	Red
2124	28 in. (71cm)	Caribbean Blue
2125	28 in. (71cm)	Raspberry
2126	28 in. (71cm)	Lilac

The Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

For more information, call the 3M Health Care Helpline at **1-800-228-3957**, or visit us at [www.littmann.com](http://www.littmann.com)

For healthcare professional use only.  
Registered under Act 737.  
MDA Reg. No. : GMD69864504618A and GMD65331541218A



#### Health Care

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
U.S.A.  
1 800 228-3957

[www.littmann.com](http://www.littmann.com)

#### 3M Canada Inc.

Post Office Box 5757  
London, Ontario N6A4T1  
Canada  
1 800-563-2921

Outside of USA, please contact your 3M Representative



*40% Pre-consumer waste paper  
10% Post-consumer waste paper*

Printed in U.S.A.  
Copyright © 3M 2005  
All Rights Reserved  
70-2008-3503-4