


3M™ Littmann® Master Cardiology Stethoscope



3M Medical
from the 3M Health Care family



3M™ Littmann® Master Cardiology Stethoscope



A tunable diaphragm and superb craftsmanship offer the ultimate performance in a superior single-sided stethoscope. The Littmann Master Cardiology stethoscope meets the exacting standards of health care professionals who demand the highest acoustic response.

Features

- Tunable diaphragm conveniently alternates between low- and high-frequency sounds without the need to turn over the chestpiece
- Handcrafted solid stainless steel chestpiece
- Two-in-one tube design eliminates tube rubbing noise
- Portable convenience of a single-sided chestpiece
- Excellent acoustic seal and comfortable fit with patented 3M™ Littmann® Snap Tight Soft-Sealing eartips
- Comfortably angled, anatomically correct headset
- Patient-friendly, nonchill rim and diaphragm. Special procedures adaptor for difficult site or pediatric auscultation
- Available in a selection of classic tube colors
- Seven-year warranty

Ordering Information

Catalog #	Length	Tube Color
2159	22 in. (55,8cm)	Black
2160	27 in. (69cm)	Black
2163	27 in. (69cm)	Burgundy
2164	27 in. (69cm)	Navy Blue
2165	27 in. (69cm)	Hunter Green
2167	27 in. (69cm)	Plum
<i>Black Edition — Black-plated chestpiece and eartubes</i>		
2161	27 in. (69cm)	Black

The Littmann stethoscopes listed do not contain natural rubber latex or dry natural rubber as components in the product or packaging.

For more information, call the 3M Health Care Helpline at **1-800-228-3957**, or visit us at **www.littmann.com**

For healthcare professional use only.
Registered under Act 737.
MDA Reg. No. : GMD59426461918A



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957

www.littmann.com

3M Canada Inc.

Post Office Box 5757
London, Ontario N6A4T1
Canada
1 800-563-2921

Outside of USA, please contact your 3M Representative



40% Pre-consumer waste paper
10% Post-consumer waste paper

Printed in U.S.A.
Copyright © 3M 2005
All Rights Reserved
70-2008-3508-3