



ENTERAL
Enteral feeding pump

easymoov6
Intuitive nutrition



Value Life

A new enteral feeding pump

Designed to meet all patients needs in ambulatory and stationary medical care, the easymoov6 enteral feeding pump allows comfort and safety for patients and medical professionals in hospital and homecare.

Intuitive and available in different languages, this pump has a large screen showing all useful and relevant information.

[Intuitive]

All the information is available in one large screen: short learning curve allowing the user to manage the pump quickly.

[Quiet]

A noiseless continuous rotary peristaltic rotor allows patients to feed during the night.


[Powerful]

A continuous rotary peristaltic rotor helps delivering high viscosity food and preventing occlusion.

[Safe]

Nutrition modes as well as visual and audible alarms allow safe management of the pump by patients and/or caregivers.

An intuitive interface




easymoov6 pump			
Reference	Designation	Additional information	Qty/box
0VEPM6Xxx ⁽¹⁾	easymoov6 enteral feeding pump	Delivered with AC adapter and pole clamp	1

⁽¹⁾Xxx: version to be specified depending on the plug and language required.

A full range of giving sets

The easymoov6 pump has a full range of giving sets designed to cover all healthcare professionals' requirements and to improve safety and comfort thanks to the ENFit™ and nutrisafe2 connections.

Some of the giving sets have an anti free flow valve, allowing hydration of patients. The range includes simple and double line with the following connections:




Giving sets for easymoov6 pump					
	Reference	Proximal connection	Type of giving set	Additional information	Qty/box
Simple line giving sets	0VENEP14S	ENPlus	ENFit™	-	30
	0VENEP14SCS	Screw cap	ENFit™	-	30
	0VENEP08S	ENPlus	ENFit™	With anti reflux valve	30
	0VENEP14V	ENPlus	nutrisafe2	-	30
	0VENEP08V	ENPlus	nutrisafe2	With anti reflux valve	30
	0VENEP14VCS	Screw cap	nutrisafe2	-	30
	0VENP1P14S	500ml bag	ENFit™	-	50
	0VENP2P14S	1200ml bag	ENFit™	-	50
	0VENP1P14V	500ml bag	nutrisafe2	-	50
	0VENEP14R	ENPlus	ENFit™	With 3 way stopcock	30
	0VENEP14RCS	Screw cap	ENFit™	With 3 way stopcock	30
	0VENEP11S	ENPlus	ENFit™	No lateral connection and no christmas tree adapter	30
0VENEP11SCS	Screw cap	ENFit™	No lateral connection and no christmas tree adapter	30	
Double line giving sets	0VENEEP14S	ENPlus ; ENPlus	ENFit™	-	30
	0VENEEP14SH	ENPlus ; ENPlus	ENFit™	With anti free flow valve	30
	0VENEP2P14SHCS	1200ml bag ; screw cap	ENFit™	With anti free flow valve	50
	0VEN2P2P14SH	1200ml bag ; 1200ml bag	ENFit™	With anti free flow valve	50
Accessories	1MWPM6CHA ⁽²⁾	easymoov6 AC adapter			1
	1MWPM6ACP	easymoov6 Pole clamp			1
	0VEPM6BAG	Backpack for easymoov6 pump			1

⁽²⁾x: version to be specified depending on the plug required.

easysoft

easy**soft** is a software to help in the management of the easy**moov6** pump allowing:

- Access to the 500-events history calendar of the pump
- Customization of some pump parameters
- Help in the preventive maintenance of the pump



easy moov6 pump			
Reference	Designation	Additional information	Qty/box
1MWEASYSOFT	easy soft	Delivered as USB key with USB-Micro USB cable	1

Important: easy**soft** is available in French and English **only**.

Pump table stand

A pump table stand is now available, allowing to hang 2 nutrient bags and to stabilize the pump wherever the patient is going. The height is adjustable and once folded it fits into the easy**moov6** backpack.

easy moov6 pump			
Reference	Designation	Additional information	Qty/box
0VESUPPOMP	Pump table stand	Delivered disassembled in 3 parts	1



OBSTETRICS NEONATOLOGY ENTERAL

For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE

Reception: +33 (0)1.39.92.63.63 – Service clients France: +33 (0)1.39.92.63.81

Export customer service: +33 (0)1.39.92.64.15

Fax.: +33 (0)1.39.92.64.44 • www.vygon.com



www.vygon.com