

• Enteral Feeding System



Enteral giving set (spike)

Enteral giving set (bag)

Enteral syringe

Enteral feeding catheter



Enteral Feeding Pump EP-60

Easy to feed
Safe to use



Distributor:


Enteral Feeding System

Integrated enteral feeding system
We provide enteral giving sets,
enteral feeding catheters and
enteral syringes



Auto Anti Occlusion

Even between two meals,
feeds will not clot



4.0" Touch Screen

Touch to set parameters
Quick training and easy operation
Two user interface themes
Teaching videos inside



Robust Housing

High ingress protection level IP34
Housing protects pump
in case of feed solution leaks



Optional Heater

Heater works and stops simultaneously with the pump
Thermal insulation and safe voltage
Without additional power cord
Can reduce irritations in gastrointestinal tract



Compact Design:

Weight: 1.2 kg
Unique clamp allows
install and removal within 1 second
Anti-skid hidden handle

Superb Parameters:

Accuracy: $\leq \pm 5\%$ Delivery rate: 1-1200 ml/h
USB upload and download

Strong Battery:

Charge for 4 hours, work for 22 hours

Smart Prime:

Fast prime at 2000 ml/h
Automatically stops when prime is finished

Exclusive Design Because we care



Heater

Hidden Handle

Clamp