



HK-300 Enteral Feeding Pump




90° rotatable pole clamp

Fixation on vertical pole
or horizontal bar



Rs232 interface

2000 event logs and alarm
information can be transferred
to computer for review



Removable pump body

convenient for cleaning

Shenzhen Hawk Optical Electronic Instrument Co.,Ltd.

Address: 2B, Building No.2, Aerospace Science & Industry factory areas,
Yousong community, Longhua street, Bao'an District, Shenzhen, China.
Postal code: 518109

Tel: 0086-755-8315 1901
Email: szhk@hawkmedical.cn

Fax: 0086-755-8315 1906
Web: www.hawkmedical.cn

HK-300 Enteral Feeding Pump Specification

Specification	
Flow rate range	1 ~400 ml/h
Volume to be infused (VTBI)	0 ~9999 ml
Volume Infused (Σ)	0 ~36000 ml
Infusion accuracy	± 10%
Bolus rate	400 ml/h
Computer Interface	Rs232 (optional)
History records	2000 history records
Waterproof Level	IP×3
Mode of operation	continuous
Applicable feeding bag	Support a wide variety brands of feeding bag
Occlusion pressure	3 adjustable occlusion pressure settings: low, middle and high
Alarms	Visual and audible alarms:Door open,Occlusion,Infusion completion, infusion near over,empty, Start reminder function, Low Battery,Battery depleted,malfunction etc.
Battery	Rechargeable lithium polymer battery,7.4V,1900mAh Can operate for more than 3 hours after fully charged at 25ml/h
Battery recharge time	10h with power on, 3h with power off
Power supply	AC:100~240V,50/60Hz DC:12V ±1V
Classification	Class I, type CF
weight	≤1.4kg
Dimensions	145×100×120 mm (L×W×D)
Power consumption	25VA
Operating conditions	Temperature:5°~40°,relative humidity:10~95%, Atmosphere pressure:86~106kPa
Transport conditions	Temperature-20°~+60°,relative humidity:10~95%, Atmosphere pressure:50~106kPa
Storage conditions	Temperature-20°~+45°,relative humidity:10~95%, Atmosphere pressure:50~106kPa

Shenzhen Hawk Optical Electronic Instrument Co.,Ltd.

Address: 2B, Building No.2, Aerospace Science & Industry factory areas,
Yousong community, Longhua street, Bao'an District, Shenzhen, China.
Postal code: 518109

Tel: 0086-755-8315 1901
Email: szhk@hawkmedical.cn

Fax:0086-755-8315 1906
Web: www.hawkmedical.cn