

5LPM OXYGEN CONCENTRATOR

The first one to focus on fluctuation of an oxygen concentrator and can control the fluctuation rate within 0.1%.



Purity For Surety



ACCESSORIES

- 1 HUMIDIFIER BOTTLE
- 1 NEBULIZER KIT
- 1 CANNULA
- 1 OUTLET TUBE
- 1 DUST FILTER

Models and Functions

Model No	Max O ₂ Flow	O ₂ Purity	Volts	Watts	O ₂ Pressure	Nebulizer	Timer	Alarm	O ₂ 00.0% Display
EVOX-5S	5L/min	93%±3%	AC 220V	350W	≤70KPa	S	YES	YES	YES

Technical Specifications

System Map	Functions	EVOX-5S
Driving System	Oxygen Flow Oxygen Concentration Oxygen Output Pressure Oxygen Number Display	0.5-5L/min 93% ±3% 40-70KPa YES
Nebulizer System	Nebulizer Mode Individual Nebulizer Outlet Display Method Operating Way Nebulizing Particle Size Accessories	YES YES On LED screen By touch panel button 3.6µm±25% Nebulize nozzle and mask
Screen Display System	Display Material What Can Be Displayed ?	All by LED lights "Oxygen" light on "Nebulize" light on "Oxygen Percentage" light on "Each Running Time" light on "Accumulated Time" light on "Timer Setting Number"
Operating System	Operating Panel Oxygen Nebulizer Timer "*" and "-"	Smart Touch Start to produce Oxygen Start to run Nebulizer Setting from 0 to 3 hours (15/30 mins increase or decrease)
Safety System	Low Oxygen Alarm Oxygen Pressure Protection Alarm Pressure Faulty Alarm Power off Alarm Compressor Faulty Alarm Overheated Protection	YES YES YES YES YES YES
Electrical System	Power Consumption Auto Power off by timer Screen Sleep Mode Noise Level	350 W YES YES <46db(A)
Packing	Machine Body Size Package Size Net Weight Gross Weight Body Material	380x270x580mm 460x350x680mm 15.6kg 19 kg ABS, Glaze